

HELLO, FRESH AIR!

Hello

Three-layer
protection

skin-friendly
comfort

01.

02.

03.

Made By
Omni



OMNI

OMNI is with you Fight the virus

25g/m²

Thickened
meltblown cloth



Porous and
breathable

3layer

Triple filtering



Soft and
comfortable

GB-2010

International
standard



Low sensitivity
and non-irritating

Sorry

WE ONLY LET 300,000 PEOPLE
ENJOY OUR PRODUCTS



OMNI Intelligent Technology, IoT smart
device provider

Be careful!

Inferior masks are more terrible than pollution!

Cheap masks use unknown raw materials and have no filtering and protective effects. The harmful substances in them, on the contrary, will harm the human body.



Buying a mask is buying a meltblown cloth!

Masks that do not use meltblown cloth as
filter material are all hooligans



25g/m²

Thickened
meltblown cloth

0.5μm

Level
protection

3 levels Protect for health!

3 layer



01-

Spunbond Nonwoven

Filter large particles

02-

Meltblown nonwoven

Blocks PM2.5, dust

03-

Spunbond Nonwoven

Filter fine dust

filter
95% Filtration efficiency

SAFE

Every breath
you take needs to
be supported safely



GB 18401-2010

Formaldehyde detection
meets standards



GB 19083-2010

Microbiological and skin
irritation tests meet standards

Surrounded by softness

Non-woven inside and outside skin, as
soft as undergarments






Nonwoven
fabric

Comfort ble

Highly breathable and
non-irritating

The surface is densely packed with air
holes, comfortable and breathable, low
sensitivity and no stimulation

Three-dimensional protection



Length and width: 17.5x9.5cm

Medium size and not tighten crowded

Suitable for most face shapes, comfortable fit,
cotton elastic band, comfortable ears

We all care about you for you

De- tails

Plastic
nose bridge

Fitted
protection
design

Wide
cotton straps



Information

Model: flat

Specifications: 17.5x 9.5cm

[Product performance structure and composition]

This product consists of a mask body, a nose clip, and a mask band. The cover body is composed of three layers of filter materials; the front and back layers are polypropylene spunbond non-woven fabrics, the middle layer is polypropylene melt-blown non-woven fabric, and the nose clip is metal core plastic. The mask band is made of elastic polyester spandex. For non-sterile products,—for single use.

[Product scope and use]

It is suitable for covering the user's mouth, nose and jaw. It is used in ordinary medical environment to wear and block the mouth and nasal cavity from exhaling or spraying pollutants.

[Instructions for use]

Flatten the mask with the dark side facing outwards and the nose clip facing up. Wear it on the nose and mouth and hang the earband on the ear.

[Precautions, warnings, and informative instructions]

1. Please use within the validity period.
2. Only for one-time use, please destroy and dispose according to regulations after use.
3. The packaging is damaged and it is forbidden to use it.
- 4 Use with caution to those who are allergic to nonwovens.
- 5 Please read the instructions before use.

[Special storage conditions and methods]

Store in a relative humidity not exceeding 80%, non-corrosive gas and well-ventilated room, avoid high temperature.

