

# HELLO, FRESH AIR!

Hello

Three-layer  
protection

skin-friendly  
comfort

01.

02.

03.

Made By  
Omni

---



OMNI

# OMNI is with you Fight the virus

---

**25g/m<sup>2</sup>**

Thickened  
meltblown cloth

---



Porous and  
breathable

---

**3layer**

Triple filtering

---



Soft and  
comfortable

---

**GB-2010**

International  
standard

---



Low sensitivity  
and non-irritating

---

# Sorry

WE ONLY LET 300,000 PEOPLE  
ENJOY OUR PRODUCTS

---



OMNI Intelligent Technology, IoT smart  
device provider

**Be careful!**

# Inferior masks are more terrible than pollution!

---

Cheap masks use unknown raw materials and have no filtering and protective effects. The harmful substances in them, on the contrary, will harm the human body.



# Buying a mask is buying a meltblown cloth!

---

Masks that do not use meltblown cloth as  
filter material are all hooligans



**25g/m<sup>2</sup>**

Thickened  
meltblown cloth

**0.5μm**

Level  
protection

# 3 levels Protect for health!

## 3 layer

01-

Spunbond Nonwoven

Filter large particles

02-

Meltblown nonwoven

Blocks PM2.5, dust

03-

Spunbond Nonwoven

Filter fine dust

**95%** Filtration efficiency

# SAFE

Every breath  
you take needs to  
be supported safely



GB 18401-2010

Formaldehyde detection  
meets standards

---



GB 19083-2010

Microbiological and skin  
irritation tests meet standards



# Surrounded by softness

Non-woven inside and outside skin, as  
soft as undergarments





Nonwoven  
fabric

# Comfort ble


Highly breathable and  
non-irritating

---

The surface is densely packed with air  
holes, comfortable and breathable, low  
sensitivity and no stimulation



# Three-dimensional protection



Length and width: 17.5x9.5cm

## Medium size and not tighten crowded

---

Suitable for most face shapes, comfortable fit,  
cotton elastic band, comfortable ears

# We all care about you for you

---

## De- tails

Plastic  
nose bridge

Fitted  
protection  
design

Wide  
cotton straps



# Information

---

Model: flat

Specifications: 17.5x 9.5cm

## [Product performance structure and composition]

This product consists of a mask body, a nose clip, and a mask band. The cover body is composed of three layers of filter materials; the front and back layers are polypropylene spunbond non-woven fabrics, the middle layer is polypropylene melt-blown non-woven fabric, and the nose clip is metal core plastic. The mask band is made of elastic polyester spandex. For non-sterile products,—for single use.

## [Product scope and use]

It is suitable for covering the user's mouth, nose and jaw. It is used in ordinary medical environment to wear and block the mouth and nasal cavity from exhaling or spraying pollutants.

## [Instructions for use]

Flatten the mask with the dark side facing outwards and the nose clip facing up. Wear it on the nose and mouth and hang the earband on the ear.

## [Precautions, warnings, and informative instructions]

1. Please use within the validity period.
2. Only for one-time use, please destroy and dispose according to regulations after use.
3. The packaging is damaged and it is forbidden to use it.
- 4 Use with caution to those who are allergic to nonwovens.
- 5 Please read the instructions before use.

## [Special storage conditions and methods]

Store in a relative humidity not exceeding 80%, non-corrosive gas and well-ventilated room, avoid high temperature.

